

Forensic Reagent Colorimetric & Legal Compliance Guide

Standardized CIELAB, sRGB/CMYK profiles, Delta E tolerance margins, and NDPS commercial quantity classifications.

Reagent	Target Compound	Color Name	Swatch	HEX	CMYK Profile	ΔE Tolerance	NDPS Section
Scott Reagent	Cocaine HCl	Cobalt Blue		#1A538C	C:81 M:41 Y:0 K:45	ΔE ≤ 25.0	Sec. 21 (Commercial: 250g)
Marquis Reagent	Heroin / Morphine	Deep Violet / Purple		#541C5C	C:9 M:70 Y:0 K:64	ΔE ≤ 22.0	Sec. 21 (Commercial: 250g)
Marquis Reagent	Amphetamine / Meth	Orange to Red-Brown		#A35D2A	C:0 M:43 Y:74 K:36	ΔE ≤ 22.0	Sec. 22 (Commercial: 50g)
Duquenois-Levine	Cannabis / Resin / THC	Indigo Violet Layer		#4A346C	C:31 M:52 Y:0 K:58	ΔE ≤ 24.0	Sec. 20 (Commercial: 1kg / 20kg)
Mandelin Reagent	Methamphetamine	Dark Green		#2F6D38	C:57 M:0 Y:49 K:57	ΔE ≤ 22.0	Sec. 22 (Commercial: 50g)
Mecke Reagent	Opiates / Oxycodone	Olive to Blue-Green		#345A3E	C:42 M:0 Y:31 K:65	ΔE ≤ 22.0	Sec. 21 (Commercial: 250g)
Ehrlich Reagent	LSD / Indoles	Pink to Purple		#8D558B	C:0 M:40 Y:1 K:45	ΔE ≤ 20.0	Sec. 22 (Commercial: 0.1g)
Simon's Reagent	Methamphetamine	Bright Cobalt Blue		#236B9B	C:77 M:31 Y:0 K:39	ΔE ≤ 20.0	Sec. 22 (Commercial: 50g)